

**AED Plus® Accessories  
and Training Products**

**ZOLL®**



# Electrodes



## **CPR-D-padz®**

This revolutionary one-piece adult electrode ensures fast and accurate placement, then senses and reports the depth and rate of the chest compressions to the ZOLL AED Plus®. With a five-year shelf life, the CPR-D-padz electrode has the longest shelf life of any disposable electrode.

**Part #8900-0800-01**

## **Pedi-padz® II**

ZOLL's unique intelligent pediatric capability is built into the AED Plus. Pedi-padz II electrodes are designed for patients under the age of 8 or under 55 pounds. The AED can sense when pediatric electrodes are attached and will switch from an adult to a pediatric ECG analysis algorithm. The Pedi-padz II electrodes have a two-year shelf life.

**Part #8900-0810-01**

# Accessories

Everything you need to keep your AED Plus accessible and ready for use.



## AED Plus Black Carry Bag

Has pockets for storing an extra set of electrodes and a sleeve of batteries.

Part #8000-0802-01

## AED Plus Police Carry Bag

Black Carry Bag with "Police" on front cover.

Part #8000-0806-01



## Type 123 Lithium Batteries

(Sleeve of 10)

Part #8000-0807-01



## Large Pelican Case

Stores AED Plus, CPR-D-padz, Pedi-padz II, and sleeve of 10 batteries.

19" x 15.4" x 7.6" (48.3 cm x 39.1 cm x 19.3 cm). IP67 certified. Military Standard C-4150J.

Part #8000-0837-01



## Small Pelican Case

Stores AED Plus unit only.

16" x 13" x 6.9" (40.6 cm x 33 cm x 17.5 cm)

Part #8000-0836-01


# Training



## AED Plus Trainer2

The AED Plus Trainer2 allows students to experience what it's like to help save a life with the AED Plus. Choose from four different fixed rescue scenarios or have complete manual control of the entire rescue. A wireless remote lets you control multiple AED Plus Trainer2s. Batteries, wireless remote, and training electrodes are included.

Part #8008-0050-01



## AED Plus Demo Kit

The AED Plus Demo Kit allows you to demonstrate Real CPR Help® with a clinical AED Plus. The Demo Kit comes with a CPR-D Demo Pad and a manikin in a large carry bag. (Clinical AED Plus and simulator not included.)

Part #8000-0834-01

## AED Plus Simulator

Simulates heart rhythm for AED Plus demonstration using clinical unit.

Part #8000-0819-01



# Accessories

Everything you need to simulate a rescue with the AED Plus.



CPR-D-padz Training Electrode  
Part #8900-0804-01



Replacement Adhesive Gels for  
CPR-D-padz Training Electrode  
Part #8900-0803-01



Stat-padz® II Training Electrodes  
(Case of 6)  
For use with AED Plus Trainer2.  
Part #8900-0805-01



Pedi-padz II Training Electrodes  
(Case of 6)  
For use with AED Plus Trainer2.  
Part #8900-000861-01



AED Plus Trainer Carry Case  
Includes pocket for storing extra  
set of training electrodes.  
Part #8000-0375-01



AED Plus Trainer2 Replacement  
Wireless Remote  
Part #8008-0007



Demo Manikin  
Part #8000-0835-01



Travel Trainer  
One-piece demonstration  
tool for clinical AED Plus.  
Part #8008-0006-01



CPR-D Demo padz  
Demonstration electrode for  
use with the clinical AED Plus.  
Part #8900-5007

# Cabinets



## 9" Standard Wall Cabinet

Visibility of your AED Plus is essential so that it is easy to find in an emergency situation. The 9" alarmed wall cabinet holds the AED Plus in a carry case with an extra set of electrodes.

Part #8000-0855



## 7" Standard Wall Cabinet

The 7" alarmed wall cabinet is designed to hold the AED Plus (on a bracket without the carry case) and one spare set of electrodes.

Part #8000-0817

# Accessories

Everything you need to keep your AED Plus accessible and ready for use.



## Semi-Recessed Wall Cabinet

Alarmed AED wall cabinet holds AED Plus (on a bracket without the carry case) and one spare set of electrodes.

17.3" x 17.3" x 3.5"

44.1 cm x 44.1 cm x 8.8 cm

Part #8000-0814



## Fully Recessed Wall Cabinet

Alarmed AED wall cabinet holds AED Plus (on a bracket without the carry case) and one spare set of electrodes.

17.3" x 17.3" x 1.5"

44.1 cm x 44.1 cm x 3.9 cm

Part #8000-0811



## Mounting Bracket

Used to mount AED Plus on a wall without the carry case.

Part #8000-0809



## AED Sign Package

Includes 3D wall sign and 8.5" x 11" wall-mounted sign.

Part #8000-0825



ZOLL Medical Corporation, an Asahi Kasei Group company, develops and markets medical devices and software solutions that help advance emergency care and save lives, while increasing clinical and operational efficiencies. With products for defibrillation and monitoring, circulation and CPR feedback, data management, fluid resuscitation, and therapeutic temperature management, ZOLL provides a comprehensive set of technologies that help clinicians, EMS and fire professionals, and lay rescuers treat victims needing resuscitation and critical care. For more information, visit [www.zoll.com](http://www.zoll.com).

The Asahi Kasei Group is a diversified group of companies led by holding company Asahi Kasei Corp., with operations in the chemicals and fibers, homes and construction materials, electronics, and health care business sectors. Its health care operations include devices and systems for critical care, dialysis, therapeutic apheresis, transfusion, and manufacture of biotherapeutics, as well as pharmaceuticals, diagnostic reagents, and nutritional products. With more than 25,000 employees around the world, the Asahi Kasei Group serves customers in more than 100 countries. For more information, visit [www.asahi-kasei.co.jp/asahi/en/](http://www.asahi-kasei.co.jp/asahi/en/).

**ADVANCING  
RESUSCITATION.  
TODAY.®**

ZOLL Medical Corporation  
Worldwide Headquarters  
269 Mill Road  
Chelmsford, MA 01824  
978-421-9655  
800-348-9011  
[www.zoll.com](http://www.zoll.com)

For subsidiary addresses and fax numbers, as well as other global locations, please go to [www.zoll.com/contacts](http://www.zoll.com/contacts).